


Manual for PanGraphViewer

Yuxuan Yuan, Ricky Ma and Ting-Fung Chan*

The Chinese University of Hong Kong, Hong Kong



Version 1.0.1


2021-12-10

Table of Contents

- [Workflow for using panGraphViewer](#)
 - [Plot a genome graph](#)
 - [Check the sequence of nodes of interest](#)
 - [Show gene models that may be interrupted by some sequences/nodes in some samples](#)
- [How to use desktop-based panGraphViewer](#)
 - [Load Files for Graph Visualization](#)
 - [GFA](#)
 - [VCF](#)
 - [Selection for Plot](#)
 - [Genome Graph](#)
 - [Check node information](#)
 - [Display gene models that may be interrupted in some samples](#)
- [How to use web browser-based panGraphViewer](#)

Workflow for using panGraphViewer

Thank you for choosing panGraphViewer. Hope you can find our tool is useful for your study. In this section, we will introduce a general **workflow** (see the figure below) that can be used to explore **panGraphViewer**.



There are three key functions that we would like to emphasize.

- Plot a genome graph
- Check the sequence of nodes of interest
- Show gene models that may be interrupted by some sequences/nodes in some samples

Plot a genome graph

Depending on the purpose and preference, users can install either application on their platform.

- Before starting the application, firstly, users need to prepare an **rGFA** file, a **GFA** file or a **VCF** file that can be directly imported into the application. An output directory is also needed when processing the analyses.
- After specifying a needed file and an output directory, users can start to parse the file and then the graph information underlying would be stored in the memory. Users can adjust their purpose to decide which chromosome/sequence or a segment of it can be displayed in the display canvas by specifying the **backbone**, **chromosome/sequence** name and **coordinates**.
- Once the graph is displayed, users can use, for example, the mouse to zoom in and out to explore the graph. Hover information will show in the canvas if users move the **mouse** to the node.

Check the sequence of nodes of interest

If users find some nodes are of interest, they can check or save the sequence of the nodes.

- Usually, users can check nodes that are in variant hotspot regions or nodes falling in the gene model regions which may interrupt the genes and change the functions in some samples.
- The saved node sequences can also be imported to other tools such as **NCBI BLAST** to check the functions.


Show gene models that may be interrupted by some sequences/nodes in some samples

If users have the annotation file for the **backbone** sample, they can check gene models that are interrupted by some sequences/nodes in some samples.

- This would be of particular interest to biologists to find out if some important gene models are interrupted by some insertions, deletions, duplications or inversions in some samples.

How to use desktop-based panGraphViewer

For the **desktop-based** application, once it is opened as shown below, users can follow the following steps to explore the program.



Load Files for Graph Visualization

Here we provide two options to let users select if they want to explore a **GFA**-based graph or a **VCF**-based graph.

GFA

When selecting to plot a **GFA**-based graph, an **rGFA** file or a **GFA** file is needed. Users can click the **Select** button to specify the file and then select an output directory to perform the analyses.

Load File(s) for Graph Visualization

GFA

VCF

(Reference) Graphical Fragment Assembly (GFA) File

The GFA File

Select

Clear

Output Directory

Select

Clear

Parse GFA

Idle

Start

VCF

When selecting to plot a **VCF**-based graph, a **VCF** file is needed.

Load File(s) for Graph Visualization

GFA

VCF

The VCF File

Select

Clear

Output Directory

Select

Clear

Backbone name

Threads

4

The Backbone Fasta

Select

Clear

Parse VCF


Idle

Start


Users can click the **Select** button to specify the VCF file. An output directory is also needed to perform the analyses. Here the **Backbone name** and **Backbone fasta** are optional if

- the VCF file is a standard one with sequence/contig header and length clearly documented and users do not have a **backbone fasta**
- users have a standard VCF file but have no interest in node sequences

The program will check automatically if the given **VCF** file is a standard one once a **fasta** file is not specified. If the **VCF** file selected is not qualified, an error message will pop up.



Otherwise, a warning message will show.




By default the **Backbone name** is **backbone** and the **Threads** to parse the VCF file is **4**.

Once the files are selected and the output directory is specified, users can click the **Start** button to parse either the given **GFA** file or the **VCF** file. The program will run internally with '**Parsing... or Converting ...**' showing in the **Status** bar. Once this is completed, '**Finished in xxx s!**' will show.

Selection for Plot

After completing GFA or VCF parsing, users can go to the Selection for Plot panel to check genome graphs. Here we provide Basic settings and Sample(s) showing panels to let users customize their plots.

In the Basic settings panel, users can select a particular chromosome/sequence to display after specifying the name of the backbone sample.




The screenshot shows a web interface titled "Selection for Plot". It has two tabs: "Basic settings" (active) and "Sample(s) showing". Under "Basic settings", there are four items: "The backbone sample used in the panGraph" with a dropdown menu, "The chromosome/sequence to check" with a dropdown menu, "Start Position" with a text input field, and "End Position" with a text input field.

For the start and end positions, users can specify both of them, one of them or none. For example,

- if users specify both start and end positions, the graph will be limited to the two positions within the selected chromosome/sequence.
- if only start position is given, the program will display a graph between the start position and the end position of the selected chromosome/sequence. if
- only end position is given, the program will display a graph between 1 to the given end position within the selected chromosome/sequence.
- if both start and end positions are not given, the program will display the entire graph of the selected chromosome/sequence.

Additionally, the program will check if the given start and end positions are qualified automatically. If not, a warning message will pop up and the plot will not be generated.

For the **Sample(s) showing** panel, users can select to hide nodes from particular sample(s) if they have such a demand.



Selection for Plot

Basic settings Sample(s) showing

The remaining samples are:

- hg38
- T2T
- hg19

Sample

Remove

Add

Reset

Genome Graph

After completing the settings above, users can either click the **Plot** button in the **Plot the Graph** panel to generate the graph or customise the **Settings** first and then plot graphs.

In **Settings**, users can select **Node Shapes** to specify node shapes in either **vis.js**-based plot or **cytoscape.js**-based plot from the dropdown menus.

Node Shapes

Node shapes showing in the graph

Vis.js Graph

| | | |
|----------------|--------------|---|
| Backbone Nodes | dot | ▼ |
| SNP | dot | ▼ |
| Deletion | triangle | ▼ |
| Insertion | triangleDown | ▼ |
| Inversion | text | ▼ |
| Duplication | database | ▼ |
| Translocation | star | ▼ |

Cytoscape.js Graph


| | | |
|----------------|-----------------|---|
| Backbone Nodes | ellipse | ▼ |
| SNP | ellipse | ▼ |
| Deletion | concave-hexagon | ▼ |
| Insertion | triangle | ▼ |
| Inversion | vee | ▼ |
| Duplication | diamond | ▼ |
| Translocation | star | ▼ |

Do you want to keep the changes?

Yes

Reset

Users can also select to use either [vis.js](#) or [cytoscape.js](#) to plot graphs by the number of nodes specified in the [Graph Modification](#) panel under [Settings](#).




The image shows a 'Graph Modification' dialog box with a title bar containing a window icon, the text 'Graph Modification', and standard window controls (minimize, maximize, close). The main title is 'Graph Modification for vis.js Graph'. It contains several settings: three dropdown menus for 'Enable interaction', 'Enable layout modification', and 'Enable add/remove node', all set to 'No'; three text input fields for 'Canvas height' (980), 'Canvas width' (1400), and 'Max num. of nodes to display' (200); two more text input fields for 'Max length of seq to display (bp)' (1000000) and 'Min num. of node overlapping with a gene' (2); and a question 'Do you want to keep the changes?' at the bottom with 'Yes' and 'Reset' buttons.

| Setting | Value |
|--|---------|
| Enable interaction | No |
| Enable layout modification | No |
| Enable add/remove node | No |
| Canvas height | 980 |
| Canvas width | 1400 |
| Max num. of nodes to display | 200 |
| Max length of seq to display (bp) | 1000000 |
| Min num. of node overlapping with a gene | 2 |

By default, if users want to check graphs with ≤ 200 nodes, **vis.js** will be applied to generate graphs. Users can adjust this value (200) depending on the preference. However, if there are over **five thousands** of nodes that users want to browse in one graph, we **don't recommend** using **vis.js** as it may take a long time to load the graph file (an html file).

Depending on the screen/display size, users may also need to adjust the canvas **height and width** to make the graph fully show in the canvas if selecting **vis.js**-based plot. The hover box may run off the screen if the canvas **height and width** settings do not match with users' screen size.

Once all settings are ready, users can click the **Plot** button to generate the graphs. After the graph is shown in the display canvas, users can use the **mouse** to zoom in and out to check graph details. Users can also move the **mouse** close to a specific node to check the node information (**mouseover**, see the figure below).



From the hover box, users can capture the information, such as

- **NodeID:** the node identity
- **Resource:** the node belongs to which chromosome/sequence in which sample
- **Len:** the sequence length of the node
- **Pos:** the coordinate of the backbone node
- **Info:** if the graph is from a VCF file, it will show the type of variants and the occur position in the backbone sample. **SNP** represents single nucleotide polymorphism, **INS** represents insertion, **INV** represents inversion, **DUP** represents duplication and **TRANS** represents translocation.
- **Seq:** the first 10 nucleotides of the node sequence if it has

If users decide to use **vis.js** to plot graphs, we also provide some options in **Settings** to let users customize the plot if they don't like the default one.

- **Graph Interaction:** to enable/disable some functions when interacting with the graph

interaction

| | |
|-----------------------|-------------------------------------|
| dragNodes: | <input checked="" type="checkbox"/> |
| dragView: | <input checked="" type="checkbox"/> |
| hideEdgesOnDrag: | <input type="checkbox"/> |
| hideNodesOnDrag: | <input type="checkbox"/> |
| hover: | <input type="checkbox"/> |
| keyboard: | <input type="checkbox"/> |
| multiselect: | <input type="checkbox"/> |
| navigationButtons: | <input type="checkbox"/> |
| selectable: | <input checked="" type="checkbox"/> |
| selectConnectedEdges: | <input checked="" type="checkbox"/> |
| hoverConnectedEdges: | <input checked="" type="checkbox"/> |
| tooltipDelay: | <input type="range" value="300"/> |
| zoomView: | <input checked="" type="checkbox"/> |


- **Graph Modification:** to customise the display of the graph, for instance the **graph layout**, **edge length** and the **degree of node overlapping**.

physics

| | |
|------------------------|--|
| enabled: | <input checked="" type="checkbox"/> |
| barnesHut: | |
| gravitationalConstant: | <input type="range" value="-2000"/> |
| centralGravity: | <input type="range" value="0.3"/> |
| springLength: | <input type="range" value="95"/> |
| springConstant: | <input type="range" value="0.04"/> |
| damping: | <input type="range" value="0.09"/> |
| avoidOverlap: | <input type="range" value="0"/> |
| maxVelocity: | <input type="range" value="50"/> |
| minVelocity: | <input type="range" value="0.75"/> |
| solver: | <input type="text" value="barnesHut"/> |
| timestep: | <input type="range" value="0.5"/> |

Users can customise the graph by enabling the options and the corresponding options would be shown in the display canvas during **graph display**.


While in **cytoscape.js**-based plot, we do not provide such functions or interactions mentioned in the **vis.js**-based plot. However, one unique interaction in **cytoscape**-based plot is that users can press **Ctrl** or **command** and hold the **left click** button of the mouse to select particular node regions and then **right click** the **mouse** to show or hide nodes (see the figure below). This would be useful if users want to check a specific subgraph.



Check node information

We provide two ways to check the node information.

- One way is that users can move the **mouse** to a specific node in the display canvas and then a hover box will show the information of the node.



- Another way is that users can check the information of nodes from the **Available Nodes** panel.

Available Nodes

▼ Add Remove

List of the selected node(s)

Check

Clear

Save the Sequences

Idle

- In the dropdown menu, all nodes in the entire graph will be listed. Users can select the node that they want to check or remove it from the selected list.
- Users can also type the node id in the **List of the selected node(s)** panel by lines
- Once clicking the **Check** button, the program will show the selected nodes in **fasta** format with a header similar to the one shown in the **hover** box. However, the sequence would be the completed sequence of the node. Users can also click the **Save the Sequences** button to save the **fasta** file to the output directory. Users can import the saved sequence file to other programs, such as **NCBI-blast** for downstream analyses.

Please **note** that by default we only allow a sequence display with ≤ 1000000 bp in length. The reason to do this is that if a large sequence is needed to display, it would take a long time to show. Users can change the setting in **Settings**--> **Graph Modification** if they really want to display a large sequence.

Display gene models that may be interrupted in some samples

If users have an annotation file for the **backbone** sample, they can check gene models that may be interrupted by some sequences in some samples.

Nodes Overlap with Genes

The annotation file Select Clear


Parse the annotation file and link nodes to genes Idle Start

Node(s) overlapping with Gene: Plot

Basically, users need to

- specify the annotation file and then parse it. By default, gene models having at least 2 nodes falling in the region can be retained. Users can change this setting in **Settings** --> **Graph Modification**.
- once the parsing is completed, users can select a gene model from the dropdown menu to check nodes that fall in the selected gene region.


After clicking the **Plot** button, a plot will show in the display canvas.



- Users can enable the zoom in and out function by clicking the **Wheel Zoom(x-axis)** button on the right top panel in the canvas.
- In the canvas, users can get the selected gene model ID and the nodes falling in the gene region.
- Users can check the nodes through the **Available Nodes** panel or plot **subgraphs** using the coordinate of the gene model.

How to use web browser-based panGraphViewer

After the installation of the web-based panGraphViewer, users can follow [README](#) to open it. The login interface is like:

The image shows a web browser window displaying the 'panGraph Login' page. The page has a dark blue header with the text 'panGraph Login'. Below the header, there is a line of text: '(username: demo, password: demodemo)'. Underneath this, there are two input fields: 'Username:' followed by a text box containing a single character 'I', and 'Password:' followed by an empty text box. At the bottom of the form, there is a blue button labeled 'Log In'.

Users can use the test account ([demo](#)) or personal account to open the application.

After opening the page showing above, users can

- follow similar approaches used in the [desktop-based](#) version, users can [upload](#) and [parse](#) a GFA file, an [rGFA](#) file or a [VCF](#) file to the system plot
- pangenome graph/subgraph from the [Plot pangenome](#) panel explore the
- node information from the [Extract nodes](#) panel
- check gene models that may be interrupted by some sequences in some samples

Input

(r)GFA VCF

Uploaded (r)GFA file (*)

Select uploaded (r)GFA ▾



Parse GFA

Plot pangenome (sub)graph

Extract nodes

Plot nodes falling in a gene region

Plot parameters (*)


Select backbone ▾

Select chr ▾

Enter start (optional)

Enter end (optional)

Plot




PanGraphViewer is a tool helping to create a graph-based pangenome using a given VCF file with or without a FASTA file. It can also accept a reference graphical fragment assembly (rGFA) file from other software, such as minigraph. PanGraphViewer shows genome graph in a graphical user interface (GUI) or a webpage.


During graph display, users can customise the layout of a graph. Users can also check genes overlapping with specific genomic regions (nodes). Multiple nodes selection enable users to browse the selected nodes only or the rest nodes only.

PanGraphViewer also provides functions to check the sequence of selected nodes and enable users to save/download the corresponding sequences when right click the mouse.

PanGraphViewer provides an easy way to display a pangenome graph from rGFA or VCF files.



For the **admin** account, users can login using the steps shown in the **README** file. The interface of the login page is like:



The image shows the Django administration login page. It has a dark blue header with the text "Django administration". Below the header, there are two input fields: "Username:" and "Password:". Below the password field, there is a blue button labeled "Log in".

Once entering the login page, the admin can create accounts for either groups or users.



The image shows the Django administration site administration page. It has a dark blue header with the text "Django administration". Below the header, there is a section titled "Site administration". Under this section, there is a table with the following content:

| AUTHENTICATION AND AUTHORIZATION | |
|----------------------------------|--|
| Groups | + Add ✎ Change |
| Users | + Add ✎ Change |